

M1D4: Prepare protein expression system

02/19/2016

Want a research relationship like this?

Course 20 | Course 7 | Course 5 |
Course 10B | Course 6-7

2/25 4:30-5:30
56-614

**UNDERGRAD & GRAD
STUDENT UROP MIXER**

Present your research to a
variety of undergrads and get
a UROP!

Refreshments will be served!

RSVP at: <https://goo.gl/A69WLQ>

Email questions to
biotech-undergrad-officers@mit.edu

Sponsored by:



SPRING IS ^{ALMOST}₁ HERE!

You are invited to a BE / "Course 20"

Spring Term Open House

for faculty, staff & undergraduate students.

Friday, February 26, 2016

Between 2:30pm and 4:30pm - stop by anytime, stay as long as you like

Room 16-341

Come for refreshments, good company and BE swag!



Today in lab



- Isolate mutants DNA by Qiagen mini-prep
- Obtain BL21(DE3)pLysS in mid-log phase, make them competent, and transform with X#Z #1, X#Z #2, wt IPC, or no DNA
- Prepare mutant DNA for sequencing
- Count colonies from X#Z plate

Results & Discussion section for M1 summary

State concisely in bullet point:

Use this slide as a checklist !

1. What was the overall goal of the experiment? (introductory sentence)
2. What was your expected result?
 - What are the expected band sizes on your gel?
3. What evidence do you have that your result is correct or incorrect?
 - What controls did you perform and did they work as expected?
4. What was the result?
 - Were bands of the expected size present? Why or why not?
5. In sum, what do these data suggest or indicate? Think about how the data were used.
6. What does this motivate you to do next?


Assignments for M1



- Protein engineering summary draft
 - due by 5pm on Saturday, March 12
 - revision due by 5pm on Monday, March 28


Summary content

1. Title
2. Abstract
3. Background, Motivation
4. Figures, Results & Discussion, Interpretation
5. Implications, Future Work



- Mini presentation due on March 16
- Blog post for M1 due March 29

Since last time...



❖ Only 1 W/F team had colonies. Show me your pipetting skills!


Transformation controls & outcomes

sample	expectation / what if?	role
no DNA	no colony. What if many? ampicillin not working wrong antibiotic contamination	negative control
control NEB SDM	many. What if none/few? damaged DNA, not Amp resistant anymore bad kit killed NEB 5alpha cells	positive control for reagents, for technique
your X#Z or wt IPC	some. What if X#Z << control? damaged plasmid or low transformation toxic plasmid low [plasmid], poor SDM efficiency	experiment

Mini-prep to isolate DNA

using QIAprep kit by Qiagen

steps		contents	purpose
prepare	P1	Tris/EDTA buffer RNase	resuspend NEB 5alpha cells weaken the cell walls
lyse	P2	SDS surfactant NaOH alkaline lysis	solubilize proteins, denature DNA
neutralize	N3	acetic acid, chaotropic salt, potassium acetate	short (plasmid) renature long DNA (rest of E. Coli) precipitate
			clear lysate
concentrate		spin: bind to silica column	pellet "garbage"
wash	PB, PE	ethanol use 'Qiagen waste'	** get rid of <i>all</i> ethanol
elute		water, pH 8.0	high-purity DNA





⋮




DNA vs. protein amplification in NEB 5α vs. BL21(DE3)pLysS *E. Coli*

- “BL21” can express IPC protein
 - when induced by lactose analog...
 - ...review details on M1D5!




Transform BL21(DE3)pLysS competent cells



- can express IPC protein
- made **competent** by CaCl₂
 - Ca²⁺ ions attract both DNA and liposaccharides
 - heat shock
- in exponential growth phase
 - OD₆₀₀ = 0.4-0.8
- handle very gently, or will lyse
 - *on ice* all the time, and with chilled solutions ** **fragile**
 - not vortexed **
- Cam (chloramphenicol) resistant *E. coli* strain
 - also Amp (ampicillin) resistant if IPC insert uptaken

A few brief notes on *E. coli* growth curve

- exponential phase
 - binary fission
 - OD600 ~ 0.4 - 0.8
 - machinery ready
- OD \neq absorbance
 - **OD: optical density**
 - **scattering, turbidity**
 - **yellow cells don't absorb 600 nm light!**
 - **UV light not safe**



Today in lab

- ~~Quiz~~
- ~~Prelab~~
- Isolate mutants DNA by Qiagen mini-prep
- Obtain BL21(DE3)pLysS in mid-log phase, make them competent, and transform with X#Z #1, X#Z #2, wt IPC, or no DNA
- Prepare mutant DNA for sequencing
- Count colonies from X#Z plate



OD600




=

0.75

You have an entire week for M1D6 homework: get started now!

Tips to write Methods section

- Use subsections with descriptive titles
 - Put in logical order
 - Begin with topic sentence
- Use clear and concise full sentences
 - Avoid tables and lists
 - Passive voice expected
- Use the most flexible units
 - Write concentration rather than volume



Homework due M1D6

- Methods from M1D3-M1D5:
 - Eliminate '109 specific details
 - Report concentrations (NOT volumes)
 - Do not include details about tubes and water
 - Avoid repeating information
 - Use sub-section titles
 - Include topic sentences in each section