



# M2D7: Analyze ICP-OES data and examine yeast tolerance to metal

- Prelab
  - Review Mod2 project experiments
- Examine and compile ICP-OES data
- Perform metal tolerance test using Fet4\_mutant yeast




# Overview of Mod 2 experiments


Last lab:



**D1&D2:** Fet4 site directed mutagenesis




**D3:** W303 transformation with mutants




**D4&D5:** IF detection of Fet4 mutant expression

This lab:



**D6:** ICP-OES analysis of heavy metal uptake



**D7:** Determine tolerance of mutant W303 to metal




**D8:** Statistics and Wrap-Up


Next lab:

# Metal tolerance experiment overview

- Examine  $OD_{600}$  for your Fet4\_mutant culture
- Dilute your culture to achieve 5ml of culture at  $\sim 1.0$   $OD_{600}$
- Spike your yeast culture with 100uM Cadmium
- Incubate for 2.5 hours
- Recheck  $OD_{600}$
- Add cell suspension and BacTiter-Glo reagent to 96 well plate
- Read luminescence on platereader upstairs



# What is the BacTiter-Glo assay measuring?



- Generation of bioluminescence
- The luciferase enzyme oxidizes luciferin substrate
  - utilizes Mg<sup>2+</sup> and ATP
- Light is produced as part of the reaction


Measured **bioluminescence** has a **linear relationship with ATP**  
**ATP production** is indicative of a live, **metabolically active cell**

# ICP-OES calibration data (found in pdf)

Cd (214.439 nm)

Intensity = 7379.58566373 \* Concentration + 13.00393963

Correlation coefficient: 0.99997



Standards	Intensity	Method Concentration	Calculated Concentration	% Error
Blank	9.41	0.00	0.00	N/A
Standard 1	22998.83	3.00	3.11	3.83
Standard 2	63918.04	9.00	8.66	3.78
Standard 3	205044.11	27.00	27.78	2.90
Standard 4	614049.08	81.00	83.21	2.73

- each wavelength has a calibration curve established using the known standards we generated
- Standards= 0ppm, 3ppm, 9ppm, 27ppm, 81ppm

# ICP-OES sample data (in pdf and csv file)

Sample Name: Sample 17

Date: 4/18/2023 11:22:03 AM

Rack:Tube: 2:22

Weight (g): 1

Volume (mL): 1

Dilution: 1

## Analyte Results

Label	Solution Concentration	Unit	SD	%RSD	Intensity	Calculated Concentration
Cd (214.439 nm)	12.58	ppm	0.06	0.49	92819.61	12.58 (ppm)
Cd (219.463 nm)	12.15	ppm	0.56	4.57	97.35	12.15 (ppm)
Cd (223.986 nm)	11.31	ppm	0.75	6.67	63.52	11.31 (ppm)
Cd (226.502 nm)	12.12	ppm	0.07	0.57	175365.66	12.12 (ppm)
Cd (226.742 nm)	10.62	ppm	0.70	6.59	95.29	10.62 (ppm)
Cd (228.802 nm)	12.15	ppm	0.10	0.80	52154.44	12.15 (ppm)
Cd (230.662 nm)	11.68	ppm	0.07	0.63	198.56	11.68 (ppm)
Cd (231.275 nm)	11.81	ppm	0.77	6.55	221.83	11.81 (ppm)

- Each team is a sample
- Each class has a control
  - Untransformed (WT)
  - Fet4 overexpression(Fet4)
- concentration is calculated in parts per million (ppm)
  - based on peak intensity at the listed wavelength and calibration curve

## Replicates Concentration

Label	Replicate 1	Replicate 2	Replicate 3	Units
Cd (214.439 nm)	12.51	12.64	12.58	ppm
Cd (219.463 nm)	12.77	11.99	11.69	ppm
Cd (223.986 nm)	12.13	11.14	10.65	ppm
Cd (226.502 nm)	12.12	12.05	12.19	ppm
Cd (226.742 nm)	9.89	11.29	10.68	ppm

# Mod2 project review (AKA: what are we doing again?)

- Research goal:
- How does your mutagenesis design fit into this research goal?

# Mod2 project review (AKA: what are we doing again?)

- How do your experiments fit into this goal?
  - Sequencing alignment:
  - ICP-OES:
  - BacTiter-Glo assay:



## For today:

1. Set up metal tolerance experiment
2. Examine ICP-OES data during incubation
3. Collect data on metal tolerance (last data for the module!)

## For M2D8

- Draft an outline of the Research Article discussion using the prompts on the homework section of the wiki and questions you answered for M2D7