

SUBE Annual Retreat

Saturday 9/23/23

Agenda

Le Meridien in Central Square

- Engaging workshops
- Inspiring talks
- Interactive activities
- Breakfast & lunch provided

Action Kitchen in Seaport

- Cook a chef-guided dinner as a group

Recharge and grow community!

- Connect with peers in the Biological Engineering and Biology departments
- Expand your knowledge
- Make new friends



RSVP

16 lucky participants chosen by lottery!

Enter lottery here:
<https://forms.gle/kab4TEPRNLxxpAWJ6>

<https://forms.gle/kab4TEPRNLxxpAWJ6>

M1D3: Use immunofluorescence staining to assess repair foci experiment


1. Prelab
2. Antibody staining for γ H2AX assay
3. Image coverslips

Office Hours

Jamie: 11a-12p Mondays

Noreen: 3-5pm Mondays

Becky: 10a-12p Wednesdays



Mod1 Overview

Last lab:


Treated MEF cells & fixed
with PFA

This lab:

Stain cells with fluorescent
antibodies & mount
coverslips

Next lab:

Analyze γ H2AX data & begin
CometChip




1. Use repair foci experiment to measure DNA breaks

- Examine effect of H_2O_2 +/- As on double strand DNA breaks by measuring γ H2AX foci formation





2. Use high-throughput genome damage assay to measure DNA damage

- Measure effects of H_2O_2 +/- As on DNA damage by measuring DNA migration in agarose matrix

Using immunofluorescence: γ H2AX assay to detect double-strand DNA breaks




- Histone H2AX phosphorylated at Ser139 if DSB
- Antibodies against γ H2AX (phosphorylated form)

protein of interest	 γ H2AX
primary antibody	 mouse anti- γ H2AX
secondary antibody	 goat anti-mouse
Fluorophore (conjugated to secondary antibody) exc./ em. wavelengths	 488/525 nm



Considerations for using antibodies in the lab

- Antibodies bind to specific epitopes on antigens
 - Antigens may have multiple epitopes



Primary antibodies vs secondary antibodies


- Primary antibody recognizes the antigen
 - Specific protein sequence
 - Specific conformation of protein
 - Specific state of protein (i.e. phosphorylation)
- Secondary Ab recognizes the species of the primary Ab
 - Often conjugated to tag for visualization
 - Enzyme or fluorophore
 - Amplifies signal through multiple bindings
 - Consider sample species when choosing antibodies!



Polyclonal vs. monoclonal antibodies


Polyclonal

- **How it's made:** animal (often rabbit) immunized with antigen of interest then antibodies collected from blood sera and affinity purified
- **Advantages:**
 - Less expensive and faster to produce than monoclonal
 - Multiple antibodies in one polyclonal mixture can increase antigen recognition by binding multiple epitopes
 - Especially useful for proteins with low expression
- **Disadvantages:**
 - Variability from lot to lot




Monoclonal


- **How it's made:** animal (usually mouse) immunized with antigen of interest then B cells from spleen are harvested and fused with myeloma cells to create hybridoma cell line that will continually produce single antibody clone
- **Advantages:**
 - Very consistent
 - Binds single epitope (can also be disadvantage)
- **Disadvantages:**
 - More expensive and requires animal sacrifice




Using immunofluorescence (IF): steps in protocol



Finish IF by adding DAPI, then mount slides for imaging



Mount coverslip on glass slide
with mounting media



Blue= DAPI
Green= γ H2AX

Homework discussion

Experimental Schematics

Notes on experimental schematics


Research Question: How do I make the best baguette?



Combine 1000 g All purpose flour, 1 tsp yeast, 1 tsp salt in a large mixing bowl

Add either 300g, 500g or 700g of warm water (37C)

Mix the dough in the mixing bowl, either by hand or using a rubber spatula



Turn the dough over onto a clean work surface and knead by pushing and folding the dough onto itself, creating gluten strands until smooth and elastic

Allow your dough to prove by returning it to the mixing bowl, covering with a wet tea cloth in the fridge overnight

What should be in the Title and Caption?

Title: State what is shown / represented in the schematic

Caption:


- Explain the flow of information using concise / clear language
- Expand on text shown in figure labels to eliminate excess wordiness / clutter from the figure
- Define all abbreviations / jargon / labels / symbols

Notes on experimental schematics

Research Question: How do I make the best baguette?



Combine Flour, yeast, salt & varying volumes of water into a dough and knead until smooth & elastic



Proof at 4C ON



Deflate the dough and shape into a long rectangle, let proof at RT until doubled in size



Score the dough



Bake until a deep golden brown



Collect Data via mastication and ingestion

Figure 1: Optimizing hydration for baguettes. Baguettes were made with varying levels of water. Doughs were proofed at 4C overnight, then deflated, shaped, and proved at room temperature. Quality data was collected by subject reporting on a scale of 1-5. ON = overnight. RT = room temperature.

Notes on experimental schematics...

How does Becky knit a scarf?



Buy beautiful yarn



Choose a pattern



Cast on 25 stitches



Knit



Purl

K & P & K & P


Measure size. It doesn't match the recommended gauge



Frog it!



Cast on 40 stitches



Follow the pattern until time to cast off



Block scarf to wear

Revised example:




Figure 1: Becky’s knitting process. Becky follows a specific protocol to knit a scarf. She chooses her yarn and checks the pattern before following the written pattern and blocking to complete the project. K2= knit two, P2= purl 2, YO= yarn over, K2tog= knit two together

In lab today:

1. Complete IF staining for H2AX
 1. Blocking: Noreen lecture on Research Talk
 2. Primary & Secondary Ab: View H2AX images on 7th floor microscope

HW due M1D4

Tues 9/26!

Thursday =
Comm Lab lecture, no lab

(group)

- Create an experimental schematic for the H2AX staining process

(individual)


1. Write outline for Research Talk
2. Read paper for discussion (linked on M1D4)
3. Visit Comm Lab before M1D5

Research Talk due Saturday, Sept 30

- Prepare a video of you verbally discussing your research
 - Use any device or Zoom
 - No visuals / slides
 - Do not edit / splice the video
- **Submit to Gmail account!**
 - bioeng20.109@gmail.com
 - Remember to follow file name guidelines

Presentation should be 3 min (+/- 15 sec)

- Introduce yourself
- Provide important background information
- Describe key results
 - Briefly describe critical methods used to generate important data
 - Use quantitative descriptions when discussing results
- Highlight the take-home message



What data / results should be included?

Only include information on the gamma-H2AX experiment

- How were the cells treated?
- How were the cells stained?
- How were the data analyzed?
- What are the results?

Review assignment description on wiki

Category	Elements of a strong presentation	Weight
Introduction	<ul style="list-style-type: none">• Introduce yourself and the research• Summarize the background information necessary to understand the research• State the research question	25%
Methods & Data	<ul style="list-style-type: none">• Provide ONLY the method information necessary to understand the results• Give complete and concise explanations of the results• Relate the results to the central question	25%
Summary & Conclusions	<ul style="list-style-type: none">• Highlight the key finding(s) relevant to the central question / hypothesis	25%
Organization	<ul style="list-style-type: none">• Give a logical, easy-to-follow narrative• Include transition statements	15%
Delivery	<ul style="list-style-type: none">• Show confidence / enthusiasm and speak clearly• Use appropriate language (technical or informal, as appropriate)• Be mindful of the time limit (3 minutes +/- 15 seconds!)	10%

The Research talk will be graded by Dr. Noreen Lyell with input from Dr. Becky Meyer and Jamie Zhan.