


- **Announcements**
- **Pre-lab Lecture**
  - ❖ Plasmid Overview
  - ❖ Restriction Enzymes Intro
  - ❖ PCR recap
  - ❖ Safety + Technical Tips
- **Lab Practical (~45 min)**

# Announcements

- BE seminar series:
  - Thursdays at 4:05 pm in 32-141
  - First seminar is Sept. 24<sup>th</sup>
  - Full schedule linked from BE website
- Introducing... Michelle, your TA for Module 1

# Plasmid Overview: pCX-EGFP



— extrachromosomal,  
ds circular DNA  
why? gene expression in a cell

■ EGFP → ORF, CDS

create  $\Delta 5'EGFP$

(later, combine w/  $\Delta 3'EGFP$  → recombine)

- ORI - origin of replication  
(2 on this plasmid)

near ORF → promoter / polyA  
introns / enhancer

# Intro to Restriction Enzymes



cut w/ EcoRI



EcoRI endonuclease


↓  
cuts up DNA

palindromic

created sticky ends

next time: how used  
in cloning

# Designing PCR primers




coding strand

primers: → ←

flap: useful for adding  
rest. sites

# PCR Process



MAE  
←

Melt → 95°C

Anneal → 50's-60's

Extend → 72°C

depends on TM of primer

\* ignore flap, truncation

↳ try on your own

→

- too long

- template for desired product

# PCR<sub>reaction</sub>

\* no template control → spurious products, contamination

Component	Function
DNA polymerase (Taq)	catalyzes DNA elongation
dNTPs	material for elongation
primers	selection, initiation of new DNA strands
Mg <sup>2+</sup> ; buffers (chem. environment)	provide the right
template	provide desired sequence

cofactor



# Today in Lab

- Keep PCR tubes cold!
- Write small *directly* on the PCR tubes – do not put sticky labels in the PCR machine.
- Safety and disposal for today's experiment