


M3D2 Pre-Lab:


Incubate with
ligand and apply
heat treatment for
protein
denaturation



Important Mod 3 dates

- Research proposal presentation due **Thursday, Dec 5 by 1 pm**
 - Completed in teams!
 - 12 minute presentation, submitted to Stellar
- Blog post due **Friday, Dec 6 by 10 pm**
- Mini-report due **Monday, Dec 9 by 10 pm**
 - Completed in teams!
 - 3 page word document, submitted to Stellar


Roadmap for Mod 3



Today

1. Prepare media for CETSA experiments
2. Incubate cells with ligand and harvest for CETSA

Cellular thermal shift assay (CETSA)



PROTOCOL

The cellular thermal shift assay for evaluating drug target interactions in cells

Rozbeh Jafari¹, Helena Almqvist², Hanna Axelsson², Marina Ignatushchenko¹, Thomas Lundbäck², Pär Nordlund¹ & Daniel Martinez Molina¹


¹Department of Medical Biochemistry and Biophysics, Division of Biophysics, Karolinska Institutet, Stockholm, Sweden. ²Chemical Biology Consortium Sweden, Science for Life Laboratory Stockholm, Division of Translational Medicine and Chemical Biology, Department of Medical Biochemistry and Biophysics, Karolinska Institutet, Solna, Sweden. Correspondence should be addressed to T.L. (thomas.lundback@ki.se), P.N. (par.nordlund@ki.se) or D.M.M. (daniel.martinez.molina@ki.se).

Published online 7 August 2014; doi:10.1038/nprot.2014.138

- Assesses thermal stabilization of protein in presence / absence of ligand
- As with DSF, the ΔT_m indicates protein stabilization / destabilization compared to control

CETSA cell preparation and ligand treatment

- Cells from previous laboratory session were used to seed 3 M cells for ligand treatment
- Today you will treat cells with ligand, harvest, heat, and snap freeze




Taking a closer look...

- Treat:

Ligand 14 $\begin{cases} \rightarrow 3\mu\text{M} \\ \rightarrow 30\mu\text{M} \end{cases}$] - experimental
DMSO — — Controls
rapamycin — +
- Harvest: spinning & collecting cells
- Heat: 59°C for 3 minutes
- Snap freeze: stop rxn & store


Heat causes protein denaturation

- As proteins denature, 'melted' primary structures aggregate
- Aggregates precipitate out of solution and can be removed via centrifugation



CETSA visualized via Western blot

- T_m determined by presence of protein bands following SDS polyacrylamide gel electrophoresis
- In pilot studies, 59 °C was identified as ideal temperature to assess FKBP12 in your experiment



For M3D3...

- With your laboratory partner, write a paragraph concerning the research question you would like to pursue for your research proposal
 - See prompts on wiki
 - Does not need to be your finalized research idea / project