

■ Announcements

- ❖ Discuss mid-term feedback
- ❖ FNT heads up *methods, 3Q (got rid of one)*
- ❖ Day 7: *quiz*; *staggered arrivals (~1-1.5 hrs)*
- ❖ OH: *1hr. This Mon; next Sun + Tue, 1hr; Mon, 2hr.*

■ Pre-lab Lecture

- ❖ SDS-PAGE
- ❖ Affinity purification recap
- ❖ Today in Lab (Mod 2 Day 6)

SDS-PAGE preparation

acrylamide - toxic

- You will make whole cell extracts with equal cell #s

- Based on OD_{600} reading, normalize

① $OD = 1.0$

② $OD = 0.5$



$V_{max} = 15 \mu L$

① 7.5 μL + 7.5 μL H_2O

② 15 μL

- Gel separates proteins based on size, shape, charge

make uniform

- Sample preparation

- SDS: coat proteins w/ (-) charge

in hood


- β -Me: breaks S-S bonds *tip: to in waste*

- Boiling: denature higher-order structures

- Sample Buffer has SDS, β -Me, plus glycerol, SPB dye


SDS-PAGE visualization, analysis

- Visualization: Coomassie stain (binds certain AA)
- Two ladders: visualization, quantification



Kaleidoscope

*in real-time,
pre-stained*





Unstained

respond to the stain

150 ng (@100 KDa)
750 ng (@50 KDa)
150 ng (@20 KDa)

Affinity purification

- Basis: His-tag in vector
6x, binds to metal



Today in Lab

- Lyse cell pellets in BPER
 - BSA “carrier,” protease inhibitors
 - Add lysis enzymes → 4 μL enzymed
- Run a 25 μL aliquot through SDS-PAGE
 - Two ladders also → bail these too
 - Stick with equal volumes if you have < 25
- Purify IPC protein from the rest (long!)
 - Immediately take 10 μL aliquot and measure concentration
 - The rest is stabilized w/BSA, to be titrated against calcium next time

