

M2D6: Characterizing Your Cph8 Mutant

10/30/12

Analysis of Cph8 mutant:



<u>EnvZ</u>	H243	D244	L245	R246	T247R (K+P-)	P248
wt seq	CAC	GAC	TTG	CGC	Thr = ACG	CCG
<u>Cph8</u>	H537	D538	L539	R540	T541	P542
	Kinase Dead mutant				NNY mutagenesis	
	GCC = Ala				CTY = Leu	
					CCY = Pro	
					CAY = His	
					CGY = Arg*	
					TTY = Phe	
					TCY = Ser	
					TAY = Tyr	
					TGY = Cys	
					ATY = Ile	
					ACY = Thr	
					AAY = Asn	
					AGY = Ser	again
					GTY = Val	
					GCY = Ala	
					GAY = Asp	
					GGY = Gly	
					N = G A T C	
					Y = C T	
					15 possible amino acids	
					No stops	

	U	C	A	G	
U	Phe	Ser	Tyr	Cys	U
	Phe	Ser	Tyr	Cys	C
	Leu	Ser	STOP	STOP	A
	Leu	Ser	STOP	Trp	G
C	Leu	Pro	His	Arg	U
	Leu	Pro	His	Arg	C
	Leu	Pro	Gln	Arg	A
	Leu	Pro	Gln	Arg	G
A	Ile	Thr	Asn	Ser	U
	Ile	Thr	Asn	Ser	C
	Ile	Thr	Lys	Arg	A
	Met	Thr	Lys	Arg	G
G	Val	Ala	Asp	Gly	U
	Val	Ala	Asp	Gly	C
	Val	Ala	Glu	Gly	A
	Val	Ala	Glu	Gly	G

Today:



1. Western Blot
2. Sequencing Results
3. Photograph

T/R Mutant Miller Assay

[edit]

Please enter the results of your B-gal assay and your mutation (after sequencing)

Team Color	Mutant 1 B-Gal Average (Dark)	Mutant 1 B-Gal Average (Light)	Mutant 1 T541X; X =	Mutant 2 B-Gal Average (Dark)	Mutant 2 B-Gal Average (Light)	Mutant 2 T541X; X =	Wild Type System B-Gal (Dark)	Wild Type System B-Gal (Light)	H243A Mutant B-Gal (Dark)	H243A Mutant B-Gal (Light)
Red	1324.65	1054.65		1268.2	1106.4		944.2	330.8	71.95	66
Orange	1077.76	351.56		1228.68	458.33		1105.34	418.94	65.42	63.91
Yellow										
Green	1287.79	531.26		749.18	390.12		11331.42	3597.32	23.80	497.65
Blue	1149.27	244.80		1208.32	532.60		1573.52	537.10	91.79	80.99
Pink	1597.37	528.91		1037.31	237.98		1000.73	323.69	49.52	48.48
Purple	1208.35	1018.45		1328.9	190.95		2192.9	322.85	43.15	38.80



Western Blot analysis of Cph8 Mutant:

1. Measure OD600 of 1:10 dilution of mutants and WT - NB

2. Calculate volume for 2OD:

e.g. if OD = 0.5; need 4 mL of 1:10 or 0.4 mL of undiluted

$$2 / 0.5 = 4$$

cells

3. Lyse protein with Epicentre reagents

4. Add sample buffer and boil ←


5. Separate proteins with SDS-PAGE

6. Transfer proteins to nitrocellulose membrane for

Western Blot

Western Blot analysis of Cph8 Mutant:


cells
lysis




-β-blue
BME
-glycerol

60 min
all proteins


Protein Blot on Nitrocellulose




SDS Polyacrylamide Gel Electrophoresis



Label with Specific Antibody




Detect Antibody



Reveals Protein of Interest

Lane	Sample	Volume to load
1	"Kaleidoscope" protein molecular weight standards	10 ul
2	H6-EnvZ positive control protein	40 ul
3	wild type light sensor	40 ul
4	mutant candidate 1	40 ul
5	mutant candidate 2	40 ul
6	"Kaleidoscope" protein molecular weight standards	10 ul
7	H6-EnvZ positive control protein	40 ul
8	wild type light sensor	40 ul
9	mutant candidate 1	40 ul
10	mutant candidate 2	40 ul

2 ggs/gel




SDS-PAGE

What is the mutation? DNA Sequencing



Four dye labeled dideoxynucleotides added to each reaction

Sanger seq.

Old **New**



Add dIT



'Chain terminating reaction'

Today you will determine your Cph8 mutant sequence

What is the mutation? DNA Sequencing

Wild Type Cph8 Sequence

~1000bp

ATGGCCACCACCGTACAACTCAGCGACCAATCCCTCCGTCAGCTAG
 AAACCCTCGCCATCC..... (1200 bp)
GTGGAAATCCAAAGTGCCCTGGCCCTGAAAAAGGCGA TCGTC
AACCTCATTTTGCGCCAGGCAGAAGAATTGCATATGGCGGCTGGT
 GTTAAGCAACTGGCGGATGACCGCACGCTGCTGATGGCGGGGGTA
 AGT CACGACTTGCGCACGCCG ←
 GAGATGATGAGCGAGCAGGATGGCTATCTGGCAGAATCGATCAAT
 AAAGATATCGAAGAGTGCAACGCCATCATTGAGCAGTTTATCGAC
 TACCTGCGCACCGGGCAGGAGATGCCGATGGAAATGGCGGATCTT
 AATGCAGTACTCGGTGAGGTGATTGCTGCCGAAAGTGGCTATGAG
 CGGGAAATTGAAACCGCGCTTTACCCCGGCAGCATTGAAGTGAAA
 ATGCACCCGCTGTCGATCAAACGCGCGGTGGCGAATATGGTGGTC
 AACGCCGCCCGTTATGGCAATGGCTGGGTCAAAGTCAGCAGCGGA
 ACGGAGCCGAATCGCGCCTGGTTCCAGGTGGAAGATGACGGTCCG
 GGAATTGCGCCGGAACAACGTAAGCACCTGTTCCAGCCGTTTGTC
 CGCGGCGACAGTGCGCGCACCATAGCGGCACGGGATTAGGGCTG
 GCAATTGTGCAGCGTATCGTGGATAACCATAACGGGATGCTGGAGC
 TTGGCACCAGCGAGCGGGGCGGGCTTTCCATTCGCGCCTGGCTGC
 CAGTGCCGGTAACGCGGGCGCAGGGCATGACAAAAGAAGGGTAA

primer
(~30 bp)

Area of interest

NNN

Plans for today:

1. Lab



2. Start your WB



3. Analyze sequencing data during downtime

4. If you have a mutant, pour a new photograph plate

What are you going to do next?

