M2D8: Cell viability, quantitative PCR, identification of regulatory motifs

- 1. Treat cells with DNA damaging agents and inhibitors for cell viability
- 2. Analyze qPCR results
- 3. R: Identify regulatory motifs in RNA-seq data

Extra (+ usual) Office Hours Next Week

- Tuesday April 17th 56-322 (lab):
 - 10:30am-1:30pm (Leslie)
 - 2:00pm-5:00pm (Noreen + Josephine's regular hour)
- Wednesday April 18th 56-322 (lab):
 - 10:00am-1:00pm (Josephine)
 - 2:00pm-4:00pm (Noreen) 4-5pm (Leslie)
- Thursday April 19th (56-341c), 10-11am (Josephine)
- Friday April 20th (56-341c), 4-5pm (Leslie)

Mod2 Research Report (20% of final grade)

Due Saturday 4/21 at 10pm

- Title, Abstract
- Introduction
- Methods
- Results (Figures and captions)
- Discussion
- References

Last week of Mod2!



Etoposide is a drug/chemotherapy that causes DNA double strand breaks

- Mechanism of action: forms ternary complex with DNA and topoisomerase II enzyme, prevents re-ligation of the DNA strands
- Cancer cells (quickly dividing cells) rely on topoisomerase II more than normal cells

Topo Type II = topoisomerase II enzyme



Ma, J. & Wang, M.D. Biophys Rev (2016) 8(Suppl 1): 75. https://doi.org/10.1007/s12551-016-0215-9

Measuring synthetic lethality in our parental and BRCA2-/- cell line

What is synthetic lethality?

combination of deficiencies (eg. knockout, LOF, drug) in 2 or more genes resulting in cell death. one deficiency alone does not cause death What parallel pathways are we perturbing in this experiment?

PSB HR (BRCA2-1-) NHEJ (Loperamide, Mibifradil) What is the result we will assess?

cell viability: cell titer glo-measures ATP via luminescence

Six compounds identified that target NHEJ and / or HR

		Repair activity	activity
0	Drug name	NHEJ	HR
F	Pimozide	0.28	0.55
ک ا	_operamide	0.20	0.57
<u>→</u> N	Mibefradil	0.28	0.57
E	Etoposide	0.65	0.08
S	SR 59230A	0.27	0.58
A	AMN082	0.19	0.92

- Loperamide = slows contractions of intestines, treatment for gastrointestinal ailments
- Mibefradil = blocks calcium channels, treatment for heart conditions

(Slide from Noreen's Lecture 4/5)

Goglia et al. (2014) Molecular Cancer Therapeutics 14:326-342

Synthetic lethality part 1: experiment overview

- 1. Choose miberfradil or loperamide, sign up at front bench
- 2. Induce double strand breaks (etoposide 37°C for 60min)
- 3. Remove etoposide media and incubate with appropriate concentration of miberfradil or loperamide till M2D9



qPCR (quantitative PCR) is used to detect and quantitate gene expression

- Fluorescence is a function of <u>dSDNA</u> concentration via SYBR green dye
- Initial DNA concentrations are proportional to
 <u>RNA</u> purified from cells (from which we made
 cDNA)
- We can compare <u>expression</u> of a particular gene in different conditions by measuring the abundance of the gene-specific transcript
- Expression of the gene of interest is normalized to a housekeeping gene, <u>GAPDH</u>



Calculate relative amounts of cDNA based on threshold cycle (C_T)



Calculate relative amounts of cDNA based on initial threshold cycle (C_T) **Goal**: Find relative expression of p21 to GAPDH = $\frac{[p21 \text{ mRNA}]}{[GAPDH \text{ mRNA}]} = \frac{[p21 \text{ cDNA}]_0}{[GAPDH \text{ cDNA}]_0}$ -cDNA concentration (Assuming perfect exponential amplification during each cycle of PCR: $[p21] = [p21]_0 * 2^{Cycle#}$ $[GAPDH] = [GAPDH]_0 * 2^{Cycle#}$. Fluoppicence of [dsD/A] at G_{1} , [P21] = [GAPDH] CTGAPDH $[P21] \times 2^{CTP21} = [GAPDH] \times 2^{G}GAPDH$ · @ CT, ~ fixed [ds DNA] assuming Similar lengths $\frac{[p21 cDNA]_0}{[p21 cDNA]_0} = 2^{(CTGNPDH - (TP2I))} = 2^{(CTGAPDH)}$ [GAPDH cDNA]₀

qPCR melt curve indicates the number of dsDNA products in reaction



What would cause multiple peaks? multiple products G primers created off-farget product, sequence similar > accidentally put multiple primer pairs for different genes in same well

-> Splice variants

Computational exercise—transcription factor binding site motifs

- Calculate position weight matrices
- Search public database of transcription factor binding
- Scan sequences to look for matching motifs
- Practice expectation-maximization algorithm for de novo motif discovery



Today in lab

- 1. Drug treat cells in tissue culture:
 - 1st: Pink, Purple, Platinum
 - 2nd:Red, Orange, Green, Blue
- 2. Analyze qPCR data
- 3. Complete "Transcription Factor Motifs" R exercise

HW due M2D9:

- Create figure of qPCR analysis with related RNA seq data (plot of p21 expression in various conditions), including figure title and caption
- Associated results section
- Associated discussion section