M2D1: Prepare cells for RNA purification

03/08/2018

- 1. Prelab discussion
- 2. ½ class to TC to seed cells for RNA purification
- 3. ½ group paper discussion of Dietlein et al.

Office hours

Noreen

- M 2-5pm
- in 16-317

Leslie

- W and F 4-5pm
- in 56-341c or lab

Josephine

- T 4-5pm
- R 10-11am
- in 56-341c or lab

+ Extra

Saturday, 03/10

- 10am-5pm
- in 56-302


Please email us if you can't make office hours and we will schedule a time to meet!

Sign up for journal club


- Pick 1 of 24 papers, or suggest your own
- Present M2D6 (April 3rd) or M2D7 (April 5th)
- Sign up by adding your name next to paper [LMM/TR/Color]
 - first come first serve!
 - you cannot switch paper after M2D2 (March 15th) March 20 5
 - only one T/R presenter and one W/F presenter per article

Slot	Day 6 (T/R)	Day 7 (T/R)	Day 6 (W/F)	Day 7 (W/F)
1				
2				
3				
4				
5				
6				
7				

M2: Experimental overview



Our cell lines: DLD-1 and BRCA2-/-




Note: olaparib is a PARP inhibitor (chemotherapy)

- DLD-1 = wild-type (or parental)
 - from the colon of a male with colorectal adenocarcinoma
- BRCA2-/- = mutant
 - disruption of exon 11 from BRCA2 gene
 - deficient in DNA repair (by homologous recombination)



Tissue culture sterile technique

- 70% ethanol is your BFF:
 - wipe cabinet before and after use
 - wipe everything that enters the cabinet
- Do not disturb air flow:
 - Do no block grille or slots
 - Minimize side-to-side arm movements
 - Work > 6" away from sash
 - Leave blower on
- Do not talk into incubator!
- Only open sterile items in hood



Mammalian cell culture medium

What do cells need to survive? growth/divison/viability



Defined

- RPMI 1640 (Roswell Park Memorial Institute)
 - (a lot of phosphate)
 - often used to culture lymphoid cells
 - · Salta

«Vitamins

· phenol red: pH rein Indicator





• FBS: fetal bovine serum

UN defred

- growth factors lipids cytokines cholesterd



- penicillin - streptomycin J prevent bacterial growth

Mammalian cell culture terminology

- confluence = density
 at ~80% split cells
 - splitting= Subculturing

DLD-1

Low Density ©ATCC

High Density

• seeding =

~20-40%. of a confluent culture on hew dish

General steps for splitting cells +WHY?

1. Look at cells, estimate confluence estmete growth/ viability

2. Rinse with PBS
wash debns/antitypsin agents/serum (FBS)


3. Detach cells with trypsin breaks substrate + cell adhes rons

4. Count cells Seed specific # in new vestel


5. "Seed" new culture vessel



Calculating number of cells



- Hemacytometer Oul cell Suspersion
- · Trypan blue Stains dead cells
- # cells / mL = 10,000 x
 average of 4 corners



$$\frac{16}{4} = 4 \times 10,000 \text{ cells/ml}^{=}$$

40,000 cells/ml

Today in lab:

- 1. Tissue Culture (TC)
 - 1st: Red, Orange, Yellow, & Green
 - 2nd: Blue, Pink, Purple, White, & Grey
 - ➤ Protocols printed for TC use, no need to move laptops etc.
 - > Do not wear PPE in or out of TC room
- 2. Paper discussion of Dietlein et al.
- Homework due Wednesday, M2D2
 - Sign up for journal club day (and article > wait till west week)
 - Create a single slide from Dietlein et al.
- Don't forget about Mod1 assignments!