## SUBE Annual Retreat

**Saturday 9/23/23** 

### Agenda

#### Le Meridien in Central Square

- Engaging workshops
- Inspiring talks
- Interactive activities
- Breakfast & lunch provided

## Action Kitchen in Seaport

Cook a chef-guided dinner as a group

# Recharge and grow community!

- Connect with peers in the Biological Engineering and Biology departments
- Expand your knowledge
- Make new friends





#### **RSVP**

16 lucky participants chosen by lottery!

Enter lottery here: <a href="https://forms.qle/ka">https://forms.qle/ka</a>

b4TEPRNLxxpAWJ6



M1D3: Use immunofluorescence staining to assess repair foci experiment


- 1. Prelab
- 2. Antibody staining for γH2AX assay
- 3. Image coverslips

### **Office Hours**

Jamie: 11a-12p Mondays

Noreen: 3-5pm Mondays

Becky: 10a-12p Wednesdays





### Mod1 Overview

### Last lab:

This lab:

Next lab:






- 1. Use repair foci experiment to measure DNA breaks
- Examine effect of  $H_2O_2$  +/- As on double strand DNA breaks by measuring  $\gamma$ H2AX foci formation

- 2. Use high-throughput genome damage assay to measure DNA damage
- Measure effects of H<sub>2</sub>O<sub>2</sub> +/- As on DNA damage by measuring DNA migration in agarose matrix

## Using immunofluorescence: γH2AX assay to detect double-strand DNA breaks





- Histone H2AX phosphorylated at Ser139 if DSB
- Antibodies against γH2AX (phosphorylated form)

protein of interest	A γH2AX
primary antibody	mouse anti-human anti-γH2AX
secondary antibody	★ goat anti-mouse
Fluorophore (conjugated to secondary antibody) exc./ em. wavelengths	488/525 nm

## Considerations for using antibodies in the lab


- Antibodies bind to specific epitopes on antigens
  - Antigens may have multiple epitopes





### Primary antibodies vs secondary antibodies

- Primary antibody recognizes the antigen
  - Specific protein sequence
  - Specific conformation of protein
  - Specific state of protein (i.e. phosphorylation)
- Secondary Ab recognizes the species of the primary Ab
  - Often conjugated to tag for visualization
    - Enzyme or fluorophore
  - Amplifies signal through multiple bindings
  - Consider sample species when choosing antibodies!



### Polyclonal vs. monoclonal antibodies


### **Polyclonal**

- How it's made: animal (often rabbit) immunized with antigen of interest then antibodies collected from blood sera and affinity purified
- Advantages:
  - Less expensive and faster to produce than monoclonal
  - Multiple antibodies in one polyclonal mixture can increase antigen recognition by binding multiple epitopes
    - Especially useful for proteins with low expression
- Disadvantages:
  - Variability from lot to lot


### **Monoclonal**

- How it's made: animal (usually mouse) immunized with antigen of interest then B cells from spleen are harvested and fused with myeloma cells to create hybridoma cell line that will continually produce single antibody clone
- Advantages:
  - Very consistent
  - Binds single epitope (can also be disadvantage)
- Disadvantages:
  - More expensive and requires animal sacrifice






## Using immunofluorescence (IF): steps in protocol



## Finish IF by adding DAPI, then mount slides for imaging




## Homework discussion

**Experimental Schematics** 

### Notes on experimental schematics...

How does Becky knit a scarf?




## What should be in the Title and Caption?

**Title:** State what is shown / represented in the schematic

### Caption:

- Explain the flow of information using concise / clear language
- Expand on text shown in figure labels to eliminate excess wordiness / clutter from the figure
- Define all abbreviations / jargon / labels / symbols

## Revised example:



**Figure 1: Becky's knitting process.** Becky follows a specific protocol to knit a scarf. She choses her yarn and checks the pattern before following the written pattern and blocking to complete the project. K2= knit two, P2= purl 2, YO= yarn over, K2tog= knit two together

## In lab today:

- 1. Complete IF staining for H2AX
  - 1. Blocking: Noreen lecture on Research Talk
  - 2. Primary & Secondary Ab: View H2AX images on 7th floor microscope

### HW due M1D4

Tues 9/26!

Thursday = Comm Lab lecture, no lab

### (group)

Create an experimental schematic for the H2AX staining process

### (individual)

- 1. Write outline for Research Talk
- 2. Read paper for discussion (linked on M1D4)
- 3. Visit Comm Lab before M1D5

## Research Talk due Saturday, Sept 30

- Prepare a video of you verbally discussing your research
  - Use any device or Zoom
  - No visuals / slides
  - Do not edit / splice the video

#### Submit to Gmail account!

- bioeng20.109@gmail.com
- Remember to follow file name guidelines

## Presentation should be 3 min (+/- 15 sec)

- Introduce yourself
- Provide important background information
- Describe key results
  - Briefly describe critical methods used to generate important data
  - Use quantitative descriptions when discussing results
- Highlight the take-home message



## What data / results should be included?

Only include information on the gamma-H2AX experiment

How were the cells treated?

How were the cells stained?

How were the data analyzed?

What are the results?

## Review assignment description on wiki

Category	Elements of a strong presentation	Weight
Introduction	<ul> <li>Introduce yourself and the research</li> <li>Summarize the background information necessary to understand the research</li> <li>State the research question</li> </ul>	25%
Methods & Data	<ul> <li>Provide ONLY the method information necessary to understand the results</li> <li>Give complete and concise explanations of the results</li> <li>Relate the results to the central question</li> </ul>	25%
Summary & Conclusions	Highlight the key finding(s) relevant to the central question / hypothesis	25%
Organization	<ul> <li>Give a logical, easy-to-follow narrative</li> <li>Include transition statements</li> </ul>	15%
Delivery	<ul> <li>Show confidence / enthusiasm and speak clearly</li> <li>Use appropriate language (technical or informal, as appropriate)</li> <li>Be mindful of the time limit (3 minutes +/- 15 seconds!)</li> </ul>	10%

The Research talk will be graded by Dr. Noreen Lyell with input from Dr. Becky Meyer and Jamie Zhan.